


FILEX

POSTERIOR LUMBAR MODULAR MIS FIXATION PLATFORM



FileX modular percutaneous polyaxial screws are designed with versatility in mind. Made for true percutaneous implantation, they allow for complete pedicle stabilization, and are designed with a proprietary mid-tower accordion hinge. This system is to be used in conjunction with the PICASSO II™ system.

Flexible tower provides superior usability in crowded areas

Stability ring prevents splay during implantation, reduction and set screw insertion

Laser scored and tapered housing provides a clean break every time

6.0 percutaneous **rod available in straight and curved**, and designed with a bullet tip, which allows for smooth insertion

25mm of internal threads aid in rod reduction

Dual lead polyaxial screw with an angulation of $\pm 20^\circ$



Flex Screwbody: 5.5mm Diameter

Part Number	Color	Length (mm)
126.1636	Green	35
126.1641	Green	40
126.1646	Green	45
126.1651	Green	50
126.1656	Green	55

Flex Flexible MIS Screw Tower

Part Number	Color
126.0034	Blue

Flex Screwbody: 6.5mm Diameter

Part Number	Color	Length (mm)
126.1836	Gold	35
126.1841	Gold	40
126.1846	Gold	45
126.1851	Gold	50
126.1856	Gold	55
126.1861	Gold	60

Flex Screwbody: 7.5mm Diameter

Part Number	Color	Length (mm)
126.2036	Magenta	35
126.2041	Magenta	40
126.2046	Magenta	45
126.2051	Magenta	50
126.2056	Magenta	55
126.2061	Magenta	60

Flex Screwbody: 8.5mm Diameter

Part Number	Color	Length (mm)
126.2241	Light Blue	35
126.2246	Light Blue	40
126.2251	Light Blue	45

